



ENT PRO MED[®]
Healthcare Products

more than breathing



2024 Product Catalogue

After having 18 years of experience in export of medical products, **ENTPROMED Healthcare Products Industry and Trade Inc.** was established in 2020 as the manufacturer and the international supplier of disposable medical items particularly E.N.T. post-operative nasal products with our motto **"More than breathing..."**

The headquarter and warehouse of ENTPROMED Healthcare Products Inc. are both located in Sultanbeyli region on Asian side of Istanbul – TURKIYE.

In the second half of 2020, Entpromed has got wider range of products such as MEDICAL TAPES, FIRST AID STRIPS, ELASTIC BANDAGES, GAUZE BANDAGES and COSMETIC STRIPS.

As we offer a well-organized system which assures effortless fast purchase experiences to our business partners in different countries, and due to high demand of making cooperation on various medical products from different business partners, ENTPROMED has become one of the most important disposable medical device and raw material supplier in the end of 2020.

ENTPROMED products and quality are preferred in 30 countries on 5 continents.



- Coated or Uncoated Carriers (Silk - Paper - Cotton - PE - PU - Nonwoven)
- Silicone Papers - Medical Paper - Cold Seals
- Hotmelt - Acrylic - Zinc Oxide Adhesives



Business License
by Municipality of
Istanbul



Free Sales Certificate
by MOH Türkiye



Declaration of
Conformity CE



ISO 13485:2016



ISO 9001:2015

ENT PRO MED[®]
Healthcare Products
more than breathing...

Manufacturer,
Exporter & Global Partner

Türkiye	Irak	Serbia	Libya	Lebanon
U.K.	Iran	Poland	U.A.E.	Armenia
U.S.A.	Azerbaijan	Spain	Romania	Chile
Venezuela	Greece	Vietnam	Belgium	Uruguay
Mexico	Hungary	Kazakhstan	Germany	Bolivia
Brazil	Kosovo	Turkmenistan	Ecuador	Yemen

Thermoplastic External Nasal Splints

Professional & Unique

Rigid form of Thermoplastic material has highest protection whereas the light weight is very comfortable comparing to all other conventional previously-used fixation devices.

ENTPROMED high quality Thermoplastic External Splints can be molded easily by hand within 60 seconds after soaking in hot water at 60-70C degrees.



Available in 5 colours



Description	Colour	Small	Medium	Large
Thermoplastic Kit	Beige	TH001	TH002	TH003
Thermoplastic Kit	Baby Blue	TH004	TH005	TH006
Thermoplastic Kit	Navy Blue	TH007	TH008	TH009
Thermoplastic Kit	Red	TH010	TH011	TH012
Thermoplastic Kit	Pink	TH013	TH014	TH015



1 shipping carton contains 50 packs / kits




Kit Components

- 1 piece Self-Adhesive External Nasal Splint
- 4 pieces Self-Adhesive PRE-CUT NASAL SHAPED Nonwoven Post-operative Plasters - Hypoallergenic
- 4 pieces Self-Adhesive Nonwoven Post-operative Plasters 10cm x 1,5cm - Hypoallergenic
- 1 piece Self-Adhesive Dorsal Pad
- 1 piece Alcohol Swab



Thermoplastic

For best fitting, available in 3 different sizes.



Self-Adhesive Nasal Splints

Aluminium Pro

Professional & Unique



Available in 10 colours




ENTPROMED Aluminium PRO external nasal splints can be moulded according to the shape of the nose in more detail than our normal ENTPROMED aluminium external nasal splints. Doctors switch to ENTPROMED Aluminium PRO nasal splints, especially after post-operative usage of thermoplastic nasal splints for a week. Thus, after the initial shape of the nose is fixed with the thermoplastic nose splint, the nose continues to be protected in the most detail with the best moulding by ENTPROMED Aluminium PRO nasal splint.



Description	Colour	Small	Medium	Large
Aluminium Kit	Beige	AGP002	AGP001	AGP003
Aluminium Kit	Baby Blue	AGP004	AGP005	AGP006
Aluminium Kit	Navy Blue	AGP007	AGP008	AGP009
Aluminium Kit	Pink	AGP010	AGP011	AGP012
Aluminium Kit	Red	AGP013	AGP014	AGP015
Aluminium Kit	Lilac	AGP016	AGP017	AGP018
Aluminium Kit	Green	AGP019	AGP020	AGP021
Aluminium Kit	Orange	AGP022	AGP023	AGP024
Aluminium Kit	Black	AGP025	AGP026	AGP027
Aluminium Kit	Yellow	AGP028	AGP029	AGP030

1 shipping carton contains 50 packs / kits





Kit Components

- 1 Aluminium Self-Adhesive Nasal Splint
- 1 Self-Adhesive Dorsal Pad
- 4+4 Nasal Hypoallergenic Nonwoven Self-Adhesive Plasters for Sensitive Skin
- 1 Alcohol Swab



Aluminium pro

For best fitting, available in 2 different sizes.



Small

Medium

Aluminium

External Nasal Splints

Professional & Unique

ENTPROMED Aluminium External Splints can be molded easily by hand within a few seconds. Aluminium structure has the highest protection whereas the light weight is very comfortable comparing to all other conventional previously-used fixation devices.



Available in 10 colours




Description	Colour	Small	Medium	Large
Aluminium Kit	Beige	AG002	AG001	AG003
Aluminium Kit	Baby Blue	AG004	AG005	AG006
Aluminium Kit	Navy Blue	AG007	AG008	AG009
Aluminium Kit	Pink	AG010	AG011	AG012
Aluminium Kit	Red	AG013	AG014	AG015
Aluminium Kit	Lilac	AG016	AG017	AG018
Aluminium Kit	Green	AG019	AG020	AG021
Aluminium Kit	Orange	AG022	AG023	AG024
Aluminium Kit	Black	AG025	AG026	AG027
Aluminium Kit	Yellow	AG028	AG029	AG030

1 shipping carton contains 50 packs / kits



Kit Components

- 1 piece Self-Adhesive External Nasal Splint
- 4 pieces Self-Adhesive PRE-CUT NASAL Splint
- 4 pieces Nonwoven Post-operative Plasters - Hypoallergenic
- 4 pieces Self-adhesive Nonwoven Post-operative Plasters 10cm x 1,5cm - Hypoallergenic
- 1 piece Self-Adhesive Dorsal Pad
- 1 piece Alcohol Swab



Aluminium

For best fitting, available in 3 different sizes.



ENT PROMED[®]
Healthcare Products
more than breathing...


Aluminium STINGRAY External Nasal Splints

Professional & Unique

ENTPROMED Aluminum External Splints can be molded easily by hand within a few seconds. Aluminium structure has the highest protection whereas the light weight is very comfortable comparing to all other conventional previously-used fixation devices.



Splint



Alcohol Swap




Dorsal Ped



59 mm



Aluminium GLITTER External Nasal Splint



Available in 5 colours in 3 sizes



Internal Nasal Splints Silicone with airway

- Airways allow post-operative nasal breathing
- Silicone material prevents adhesion on skin
- Suture holes
- Made of Soft Silicone for easy insertion and removal

Indications

- Septoplasties
- Sinusotomy
- Rhinoplasties
- Turbinectomies
- Allergic rhinitis surgeries
- Nasal valve collapse surgeries



Internal Nasal Splints

Silicone without airway

- Silicone material prevents adhesion on skin
- Suture holes
- Made of Soft Silicone for easy insertion and removal



Indications

- Septoplasties
- Sinusotomy
- Rhinoplasties
- Turbinectomies
- Epistaxe surgeries
- Hypophyisic surgeries

PVA

Nasal Dressings

Professional & Unique

Indications

Epistaxis and post-operative hemostasis of nasal cavity and paranasal sinus.

Product Properties

Made of PVA (Polyvinylalcohol) expansible sponge. Soft, comfortable and conforming use. Vacuum processed sterile packaging ensuring long and safe protection. Easily trimmed with scissors or by hand.

Micro pore size minimizes tissue ingrowth. Capable to absorb 25 times of its own weight. Expansion capability under mucosal pressure. Biocompatible.



* With string

* Without string

ENT PROMED



ENT PROMED

Made in Turkey

Sterile R

- Single use - Gamma Sterilization
- Do not use if the package is damaged
- Store in a cool and dry place between 0C and 40C

ENT PROMED SAĞLIK ÜRÜNLERİ SAN. VE TİC. A.Ş.
Add: Mah. Çarşamba Sok. No: 44. Subarbaylı - İSTANBUL - TÜRKİYE
info@entpromed.com - Tel: +90 545 398 3010 - www.entpromed.com

PVA Nasal Dressing
RSF : EPN005
10cm x 1,5cm x 2,5cm
With String



8 684000 857325

LOT : 2312201

EXP : 11 / 2028



ENT PRO MED[®]
Healthcare Products
more than breathing...

The Nostril Retainer is intended for the prophylaxis of stenoses and scarring contractures. Single-use.

Indications :

- Preoperative Use: “Nostril Retainer” can be applied in the correction of alar cartilage in the cleft lip nose.
- Post-operative Use: “Nostril Retainer” is used to shape, correct and restructure the nostrils after cleft lip / palate or rhinoplasty and nose enlargement surgeries.

Nostril Retainer

Professional & Unique

Nostril Retainer is made of soft silicone that does not harm the human body. It is used to correct or reconstruct the shape of the nose by placing it inward through the nostrils. The raw material of the device is “silicone gum”. This product prevents the deformation or collapse of the nostrils after the nose aesthetic



surgery, nasal valve operation and for cleft lip or palate surgery.

It enables nasal breathing of the patient in immediate postoperative period, not obstructing the air flow. It is easy to clean and is designed for continuous use, for a period of approximately 6 months, and may be replaced by another nasal modeler of a larger size, according to the growth of the nose.¹³



For best fitting, available in 3 different sizes.

www.entpromed.com



ENT PRO MED[®]

Healthcare Products

more than bread

External Nasal Sling

Unique Patent Article
Absorbent Nose Bandage
For All Nasal Bleedings



Description	Colour	Standart Size
Sling Bandage	White	NSB02



ENT PROMED[®]
Healthcare Products
more than breathing...



catalogue



web



whatsapp

ENTPROMED SAĞLIK ÜRÜNLER SAN. VE TİC. A.Ş.

Adil Mah. Ekol Cad. Cesaret Sok. No:4A, 34935 Sultanbeyli, İstanbul – TÜRKİYE

Tel : +90 216 592 3010, info@entpromed.com

www.entpromed.com