





NETWORK

MEDICAL PRODUCTS



NETWORK

ENT PRODUCTS

ENT Single Use Disposable Products

NETCELL

ENT PVA Sponge Products



NETCELL NASAL PACK

- For use after septal, turbinate and rhinoplastic surgery.
- Provides gentle, equal pressure within the nasal cavity when absorbing post-operative bleeding.
- Available with an airway tube which will enhance patient comfort.
- Slimpack developed for use in conjunction with a nasal splint following septal surgery.
- Layered slimpack permits the user to customise the thickness of the pack.



NETCELL NASAL PACK WITH TUBE

The Netcell Nasal Packing can be supplied with an integral airway tube. The airway tube benefits the patient through easier breathing and improved patient comfort, it improves venting and reduces the risk of obstruction.

Supplied Sterile

Product	Dimensions L x Th x H (cm)	Pack	Code
Nasal Pack 3.0 cm	3.0 x 1.5 x 2.0 cm	Box of 20	10-095
Nasal Pack 4.5 cm	4.5 x 1.5 x 2.0 cm	Box of 20	10-100
Nasal Pack with Tube 4.5 cm	4.5 x 1.5 x 2.0 cm	Box of 20	10-100T
Nasal Slimpack 4.5 cm	4.5 x 1.0 x 2.0 cm	Box of 20	10-101
Nasal Pack 6.0 cm	6.0 x 1.5 x 2.0 cm	Box of 20	10-105
Nasal Pack with Tube 6.0 cm	6.0 x 1.5 x 2.0 cm	Box of 20	10-105T
Nasal Slimpack 6.0 cm	6.0 x 1.0 x 2.0 cm	Box of 20	10-106
Nasal Pack 8.0 cm	8.0 x 1.5 x 2.0 cm	Box of 20	10-110
Nasal Pack with Tube 8.0 cm	8.0 x 1.5 x 2.0 cm	Box of 20	10-110T
Nasal Slimpack 8.0 cm	8.0 x 1.0 x 2.0 cm	Box of 20	10-111
Nasal Slimpack layered 8.0cm	8.0 x 1.0 x 2.0 cm	Box of 20	10-111L
Nasal Slimpack with Tube 8.0 cm	8.0 x 1.0 x 2.0 cm	Box of 20	10-111T
Nasal Pack 10.0 cm	10.0 x 1.5 x 2.0 cm	Box of 20	10-115
Nasal Slimpack 10.0cm	10.0 x 1.0 x 2.0 cm	Box of 20	10-116



NETCELL SINUS PACK

- Developed as a post-operative packing for F.E.S.S.
- Designed to prevent lateralization of the middle turbinates post-operatively.
- Can be used to pack the Osteomeatal complex to leave the nasal passage clear for easier breathing.

Supplied Sterile

Product	Dimensions L x Th x H	Pack	Code
Sinus Pack Standard 2.5 cm	2.5 x 1.20 x 2.0 cm	Pack of 2/box of 10 Packs	20-200
Sinus Pack Standard 3.5 cm	3.5 x 1.20 x 1.2 cm	Pack of 2/box of 10 Packs	20-210
Sinus Pack Medium 3.5 cm	3.5 x 0.90 x 1.2 cm	Pack of 2/box of 10 Packs	20-211
Sinus Pack Small 3.5 cm	3.5 x 0.60 x 1.2 cm	Pack of 2/box of 10 Packs	20-212



NETCELL EAR-WICK

- Ideal for delivering antibiotics when treating Otitis Externa.
- Instant patient relief when hydrated with medicated otic solution.
- Available as fenestrated to enable aeration

Supplied Sterile

Product	Dimensions	Pack	Code
Ear Wick 15 mm	9 mm dia x 15 mm	1 per pack / box of 50 packs	30-300
Ear Wick 12 mm	7 mm dia x 12 mm	1 per pack / box of 10 packs	30-301
Ear Wick Fenestrated 15 mm	9 mm dia x 15 mm	1 per pack / box of 50 packs	30-305
Ear Wick 24 mm	9 mm dia x 24 mm	1 per pack / box of 10 packs	30-310
Ear Wick Fenestrated 24 mm	9 mm dia x 24 mm	1 per pack / box of 10 packs	30-315



NETCELL POST-OP EAR PACK

- Designed to prevent stenosis of the canal following surgery.
- Will not expand lengthwise and is suitable to support a graft following tympanoplasty.

Supplied Sterile

Product	Dimensions	Pack	Code
Ear Pack (Strung) small 15 mm	9 mm dia x 15 mm	1 per pack / box of 10 packs	30-320
Ear Pack without string 24 mm	12 mm dia x 24 mm	1 per pack / box of 10 packs	30-340



EPISTAXIS BALLOON CATHETER

- The most effective control of posterior epistaxis.
- Enhanced control of Intra Nasal Bleeding with the dual balloon system.
- Easy placement assisted by the bevelled posterior rim design.
- Airway management and patient breathing much assisted through the integral airway tube.
- Patient comfort is helped by using soft, medical grade silicone.
- Inflatable balloon cuffs conform to nasal cavity.
- Balloon size: 30cc and 10cc

Supplied Sterile

Description	Length	Pack	Code
Epistaxis Balloon Catheter Large	120 mm	Box of 1	72-3000
Epistaxis Balloon Catheter Small	97 mm	Box of 1	72-3001

NETWORK NASAL SPLINTS


- Internal and external nasal splints.
- Airway splints. Provides septal support, airway management and easier breathing.



EXTERNAL NASAL SPLINT: MALLEABLE

Supplied Non-Sterile

Description	Pack	Code
External Nasal Splint, Malleable, Small	10	72-3020
External Nasal Splint, Malleable, Medium	10	72-3021
External Nasal Splint, Malleable, Large	10	72-3022



EXTERNAL NASAL SPLINT: THERMOPLASTIC

Supplied Non-Sterile

Description	Pack	Code
External Nasal Splint, Thermoplastic, Small	10	72-3025
External Nasal Splint, Thermoplastic, Medium	10	72-3026
External Nasal Splint, Thermoplastic – Custom Cut, Large	10	72-3027



AIRWAY SPLINT: SILICONE

Supplied Sterile

Description	Pack	Code
Airway Splint – Pair left and right	10 (Pairs)	72-3030



BIVALVE SPLINT: FLUOROPLASTIC

Supplied Sterile

Description	Pack	Code
Bivalve Splint, Standard Thin	10	72-3032
Bivalve Splint, Standard Thick	10	72-3033
Bivalve Splint, Large Thin	10	72-3034
Bivalve Splint, Large Thick	10	72-3035



SEPTAL BUTTON: SILICONE


Supplied Sterile

Description	Pack	Code
Septal Button, Diameter 32mm	Single	72-3040







VENTILATION TUBES

Dimensions key: A = Internal Diameter B = Inner Flange Diameter



Supplied Sterile

Description	Material	A x B x C	Pack	Code
 T-Tube (12mm)	Silicone	1.1 x 9.0 x 12.0 mm	10 per box	74-100
 T-Tube (9mm)	Silicone	1.1 x 9.0 x 9.0 mm	10 per box	74-101
 T-Tube (6mm)	Silicone	1.1 x 9.0 x 6.0 mm	10 per box	74-102
 Paparella Tube II	Silicone	1.5 x 4.5 x 1.1 mm	10 per box	74-103
 Paparella Tube I	Silicone	1.0 x 2.1 x 1.1 mm	10 per box	74-104
 Shepard Grommets with wire	Fluoroplastic	1.1 x 2.3 x 1.5 mm	10 per box	74-106
 Shepard Grommets with wire	Fluoroplastic	1.1 x 2.3 x 1.5 mm	2 per pouch; 20 pouches	74-107
 Donaldson	Silicone	1.1 x 2.4 x 0.8 mm	10 per box	74-108
 Donaldson with tab	Silicone	1.1 x 2.4 x 0.8 mm	10 per box	74-109
 Shepard Grommets without wire	Fluoroplastic	1.1 x 2.3 x 1.5 mm	10 per box	74-110
 Shah without wire	Fluoroplastic	1.1 x 3.1 x 1.5 mm	10 per box	74-111
 Shah with wire	Fluoroplastic	1.1 x 3.1 x 1.5 mm	10 per box	74-112
 Armstrong	Fluoroplastic	1.1 x 2.6 x 4.0 mm	10 per box	74-115
 Shepard with tail	Fluoroplastic	1.1 x 2.4 x 1.5 mm	10 per box	74-117
 Donaldson	Fluoroplastic	1.1 x 2.3 x 0.8 mm	10 per box	74-120
 Reuter Bobbin	Fluoroplastic	1.0 x 2.5 x 0.9 mm	10 per box	74-121
 Collar Button	Fluoroplastic	1.2 x 2.9 x 1.3 mm	10 per box	74-122
 Collar Button	Fluoroplastic	1.0 x 2.5 x 1.2 mm	10 per box	74-123